

TriCore offers molecular assays for non-CML and non-PV myeloproliferative neoplasms. These tests may be ordered as a reflex along with *JAK2* V617F mutation analysis or individually. If ordered with the *JAK2* V617F, only one EDTA drawn tube is required.

TEST NAME	CPT CODE	METHODOLOGY	TESTING PERFORMED	REPORTED
Calreticulin (CALR) Mutation Analysis, Myeloproliferative Neoplasm	81479, G0452-26	FRAGMENT ANALYSIS	TUESDAY	Within 7 days
MPL Mutation Analysis	81402, G0452-26	PYROSEQUENCING	TUESDAY	Within 7 days

As a valued customer, the inclusion of these tests in our menu provide:

- A faster turnaround time for all 3 tests
- Single integrated report that will include results and interpretation
- One testing facility where all results can be obtained with submission of one specimen

Please use the Hematopathology requisition to order these tests located on tricore.org. Refer to our Test Directory for collection, processing and shipping information.

Should you have any questions about these tests, please contact:

Dr. Mohammad Vasef at 505-938-8468 or 1-800-245-3296 ext. 8468

Dr. David Czuchlewski at 505-938-8467 or 1-800-245-3296 ext. 8467

COLLECTION									
Specimen	Whole Blood or Bone Marrow Aspirate								
Collection Containers	Lavender-top tube (EDTA)								
Collection Amounts	Preferred 5.0 mL Minimum 2.0 mL								
PROCESSING									
Specimen Processing Instructions	CRITICAL: Do not open or share tube.								
Rejection Criteria	Frozen samples and paraffin-embedded tissues are unacceptable.								
SHIPPING/TRANSPORT AND STORAGE									
Stabilities/Storage: (Collection to initiation of testing)	<table> <tr> <td>Temperature</td> <td>Stability</td> </tr> <tr> <td>Ambient</td> <td>5 day(s)</td> </tr> <tr> <td>Refrigerated</td> <td>5 day(s)</td> </tr> <tr> <td>Frozen</td> <td>Unacceptable</td> </tr> </table>	Temperature	Stability	Ambient	5 day(s)	Refrigerated	5 day(s)	Frozen	Unacceptable
Temperature	Stability								
Ambient	5 day(s)								
Refrigerated	5 day(s)								
Frozen	Unacceptable								
Shipping Instructions	Ship ambient and overnight for Monday-Saturday receipt								