

**Effective immediately**, a new umbilical cord drug screen, UMBQL is available at TriCore. (UCDS12 was discontinued 6/26/19.) The new screen has the following benefits:

- UMBQL includes buprenorphine and fentanyl
- decreased turnaround time
- chain of custody form is no longer needed

Tests:

<b>UMBETG</b> UMBILICAL ETHYL GLUCURONIDE	<b>UMBQL</b> UMBILICAL DRUG PANEL QUALITATIVE	<b>UMBTHC</b> UMBILICAL THC
ethyl glucuronide	amphetamines barbiturates benzodiazepines buprenorphine cannabinoids cathinones cocaine fentanyl meperidine methadone opiates oxycodone phencyclidine propoxyphene tramadol	tetrahydrocannabinol

Hospital ordering codes:

<b>LOVELACE HOSPITALS</b>	<b>PRESBYTERIAN HOSPITALS</b>	<b>UNMH</b>
UMBETG	Umbilical Cord Drug Screen Panel Qualitative	UMBETG
UMBQL		UMBQL
UMBTHC		UMBTHC

## ETHYL GLU UMBIL, QL

Test Code.....UMBETG	Testing Performed.....Tuesday, Thursday, Saturday
CPT Code.....80321	Reported.....Within 2-5 days
Methodology.....Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry	

### COLLECTION

Collection Instructions	Caution must be used when collecting specimen, to ensure no ethanol-containing personal care products (i.e., hand sanitizers, wipes, mouthwash) are used directly on the specimen or nearby during collection.  Drain and discard any blood. Rinse the exterior of the cord segment with normal saline or water. Pat the cord dry and transport at least 8 inches of umbilical cord in a routine urine collection cup.
-------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## THC UMBIL QUAL

Test Code.....UMBTHC	Testing Performed.....Daily
CPT Code.....80349	Reported.....Within 2-4 days
Methodology.....Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry	

### COLLECTION

Collection Instructions	Drain and discard any blood. Rinse the exterior of the cord segment with normal saline or water. Pat the cord dry and transport at least 8 inches of umbilical cord in a routine urine collection cup.
-------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## UMBIL DRUG PNL QUAL

Test Code.....UMBQL	Testing Performed.....Daily
CPT Code.....80307	Reported.....Within 2-4 days
Methodology.....Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry	

### COLLECTION

Collection Instructions	Drain and discard any blood. Rinse the exterior of the cord segment with normal saline or water. Pat the cord dry and transport at least 8 inches of umbilical cord in a routine urine collection cup.
-------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### THE FOLLOWING INSTRUCTIONS APPLY TO THE THREE TESTS ABOVE

Specimen	Tissue (umbilical cord section)
Collection Containers	Urine container, sterile
Collection Amounts	At least 8 inches, approximately the width of a sheet of paper
Rejection Criteria	Cords soaking in blood or other fluid. Formalin fixed. Tissue that is obviously decomposed.

### SHIPPING/TRANSPORT AND STORAGE

Stabilities/Storage (Collection to initiation of testing)	<b>Temperature</b>	<b>Stability</b>
	Ambient	1 week
	Refrigerated	3 weeks
Storage Instructions	Store refrigerated	
Shipping Instructions	Ship refrigerated	