




Teal outline indicates compounds that are confirmed at TriCore Reference Laboratories with test code MBENZ.  
Times are elimination half-lives.



Screen cut-off: 200 ng/mL

LC-MS/MS Confirmation cut-off: 10 ng/mL for all compounds

### Notes

- The UBENZ test code screen is not sensitive to 7-Aminoclonazepam, and does not detect metabolic glucuronide forms of Lorazepam or Oxazepam.
- Therapeutic doses of oxaprozin (DAYPRO) may cause a false positive result.